# Luxatemp<sup>®</sup>-Glaze & Bond

CE 0482

091143/10.2007

## Gebrauchsinformation

Deutsch

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## Instructions for use

## English

Luxatemp-Glaze & Bond is a light-cure one-component material based on multifunctional methacrylates for sealing the surface of and repairing temporary crowns and bridges. The varnished surface is high-gloss, protects against discolouration and produces a highly natural appearance.

### Indication

Surface sealing and repair of temporary crowns and bridges

### Recommended use

#### Surface varnishing

1. Apply a thin layer of the material to the corresponding surfaces without air bubbles.
2. Cure each section for 20 seconds with a halogen lamp or cure the entire restoration for at least 90 seconds in the light oven.

**To ensure that no oxygen inhibited layer forms, the material should be light cured within 2:00 minutes after use.**

#### Priming for repair

1. Roughen the surface to be repaired with a bur.
2. Apply a thin layer of the material to the corresponding surfaces without air bubbles.
3. Expose to a halogen lamp for 20 seconds
4. Apply freshly mixed temporary material, e.g. Luxatemp, or a light-cure composite, e.g. LuxaFlow, to the corresponding surfaces. Follow the instructions for use of the material employed.

### Please note

- Light curing units should have an output of 450 nm and should be checked regularly. Place the light as close as possible to the material.
- Luxatemp-Glaze & Bond bonds to teeth and sheets made of polycarbonate or polyacetate.
- Luxatemp-Glaze & Bond cannot be used as a dentine or enamel adhesive.
- Avoid contact with materials containing eugenol.
- Close the bottle immediately after every use.

### Contraindication

Do not use in case of known allergy against one of the components (especially methyl methacrylates).

### Storage

Store in a dry place below 25 °C.

Do not use after the expiry date.

### Composition

Multifunctional acrylates, methyl methacrylates, catalysts, stabilizers, additives

### Presentation

#### Pack

1 bottle with 5 ml liquid, 25 brushes

REF 211927

#### Accessories

50 brushes

REF 110041

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

Keep away from children. For dental use only.

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